



# TECHNICAL DATA SHEET

**Product Name: DISINFECTANT WIPES**

**Product Number: MARC 350**

## DESCRIPTION

**MARC 350 Disinfectant Wipes** is a non-woven disposable cloth containing a stable, low pH formulated disinfectant and deodorant for use on hard, non-porous surfaces where control of cross-contamination is required. Note: May be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection. Not to be used on food contact surfaces. Effective against the following pathogens in 2 minutes at 20° C/68°F on hard, non-porous surfaces when used as directed: Staphylococcus aureus, Pseudomonas aeruginosa, Salmonella chlorosis, Escherichia coli 0157:H7, Herpes simplex virus 2/genital herpes virus, HIV-1/AIDS virus (at an exposure time of 1 minute), Trichophyton mentagrophytes/Athlete's foot fungus (at an exposure time of 5 minutes), Myco-bacterium tuberculosis, Methicillin Resistant Staphylococcus aureus (MRSA), Vancomycin Resistant Enterococcus faecalis (VRE), Influenza A. See label for complete HIV (AIDS virus) information, and other information, uses and precautions.

## APPLICATIONS

Areas of use include hospitals, intensive care units, surgery, recovery, anesthesia, X-ray, Cat Lab., orthopedics, newborn nursery, respiratory therapy, emergency medical settings, day care, schools, school buses, laboratories, clinics, nursing facilities, medical offices, dental offices, veterinary facilities, janitorial, commercial, industrial facilities, spas and health clubs, and more. For use on hard, non-porous surfaces only, such as stainless steel, Formica, glass tables, carts, baskets, counters, cabinets, telephones.

## FEATURES & BENEFITS

- Can be used anywhere on nearly any hard, nonporous surface to disinfect and deodorize.
- Safely decreases chances of cross-contamination.
- Dispenses one at a time so you only use what you need for less waste.

## SPECIFICATIONS

**FORM:** Pre-moistened Towel

**ODOR:** Citrus

**COLOR:** White Towel

**DETERGENCY:** N/A

**TOXICITY:** N/A

**WETTING ABILITY:** Excellent

**PHOSPHATES:** N/A

**STORAGE STABILITY:** 1 year+

**FREEZING POINT:** 32°F

**Size of Wipe:** 7" X 10"

**Cannister Net Weight:** 14 oz.

**FLASH POINT:** N/A

**BOILING POINT:** N/A

**SOLUBILITY IN WATER:** 100%

**SPECIFIC GRAVITY:** 1.00 - liquid

**% VOLATILE:** 100% - liquid

**EVAPORATION RATE:** N/A

**BIODEGRADABLE:** Yes

**FLAME EXTENSION:** N/A

**pH:** N/A

**Size of Cannister:** 3 ½" in diameter,  
8" high

### **ACTIVE INGREDIENTS:**

N-Alkyl (60% C14, 30% C16, 5% C12, 5% C18)

Dimethyl Benzyl Ammonium Chloride .....0.12%

N-Alkyl (68% C12, 32% C14) Dimethyl Ethyl

Benzyl Ammonium Chloride.....0.12%

Isopropyl Alcohol.....41.58%

**INERT INGREDIENTS\*\*** .....58.18%

**TOTAL**.....100.00%

\*\*Does not include weight of dry wipe.

## PACKAGING INFORMATION

50 pre-moistened towels per container, 6 containers per box.

R0621

**MID-AMERICAN RESEARCH CHEMICAL CORP.**

P. O. Box 927 • Columbus, NE 68602-0927 • Phone (402) 564-7104 • Fax (402) 563-1290 • Email: marc1@marc1.com • www.marc1.com